

## Analytical method development and validation for the estimation of Metformin and Canagliflozin using RP-HPLC

Mr. Suresh Babu, M.Pharm(Ph.D), R.Sirisha, S.Sowjanya, S.Sravani, S.Sravani

Assoc Professor Jogaiah Institute of Technology and Sciences College of Pharmacy Kallagampudi, India

\*Corresponding Author: Mr. Suresh Babu, M.Pharm (Ph.D)

Email id: Sureshbabu3377@gmail.com

### ABSTRACT

A simple and selective LC method is described for the determination of Metformin and Canagliflozin in tablet dosage forms. Chromatographic separation was achieved on a  $c_{18}$  column using mobile phase consisting of a mixture of 40 volumes of  $K_2HPO_4+NaHPO_4$  buffer, 40 volumes of acetonitrile and 20 volumes of Methanol with detection of 251 nm. Linearity was observed in the range 75-175  $\mu\text{g/ml}$  for Metformin ( $r^2 = 0.997$ ) and 24-56  $\mu\text{g/ml}$  for Canagliflozin ( $r^2 = 0.993$ ) for the amount of drugs estimated by the proposed methods was in good agreement with the label claim.

The proposed methods were validated. The accuracy of the methods was assessed by recovery studies at three different levels. Recovery experiments indicated the absence of interference from commonly encountered pharmaceutical additives. The method was found to be precise as indicated by the repeatability analysis, showing %RSD less than 2. All statistical data proves validity of the methods and can be used for routine analysis of pharmaceutical dosage form.

**Keywords:** Liquid chromatography (LC), RSD Relative standard deviation,  $R^2$  Correlation coefficient.

### INTRODUCTION

A drug includes all medicines intended for internal or external use for or in the diagnosis, treatment, mitigation or prevention of disease or disorder in human beings or animals, and manufactured exclusively in accordance with the formulae mentioned in authoritative books.<sup>1</sup>

Pharmaceutical analysis is a branch of chemistry involving a process of identification, determination, quantification, purification and separation of components in a mixture or determination of chemical structure of compounds. There are two

main types of analysis – Qualitative and Quantitative analysis.

### AIM AND PLAN OF WORK

#### Aim

To develop new RP HPLC method for the estimation of METFORMIN AND CANAGLIFLOZIN in pharmaceutical dosage form.

#### Plan of Work

- Solubility determination of METFORMIN AND CANAGLIFLOZIN various solvents and buffers.

- Determine the absorption maxima of the drug in UV–Visible region in different solvents/buffers and selecting the solvents for HPLC method development.
- Optimize the mobile phase and flow rates for proper resolution and retention times.
- Validate the developed method as per ICH guidelines.

## METHODOLOGY

### Mobile Phase

A mixture of  $K_2HPO_4+NAHPO_4$  : ACN:MEOH were prepared. The mobile phase was sonicated for 10min to remove gases and filtered through  $0.45\mu$  membrane filter for degassing of mobile phase.

### Determination of Working Wavelength ( $\lambda_{max}$ )

In estimation of drug wavelength maxima is used.. So this wavelength is used in estimation to estimate drug accurately.

### Preparation of standard stock solution of METFORMIN

10 mg of METFORMIN was weighed and transferred in to 10ml volumetric flask and dissolved in water and then make up to the mark with water and prepare  $125\mu\text{g/ml}$  of solution by diluting 1.25ml to 10ml with water.

### Preparation of standard stock solution of CANAGLIFLOZIN

5 mg of CANAGLIFLOZIN was weighed in to 10ml volumetric flask and dissolved in water and then dilute up to the mark with water and prepare  $40\mu\text{g/ml}$  of solution by diluting 0.8ml to 10ml with water.

## RESULTS AND DISCUSSIONS

### Solubility Studies

These studies are carried out at  $25^{\circ}\text{C}$

### METFORMIN

Freely soluble in methanol, water and mixed phosphate buffer.

### CANAGLIFLOZIN

Freely soluble in ethanol and methanol, and slightly soluble in acetone and very slightly soluble in water.

### Wavelength determination

In simultaneous estimation of two drugs isobestic wavelength is used. Isobestic point is the wavelength where the molar absorptivity is the same for two substances that are interconvertible. So this wavelength is used in simultaneous estimation to estimate both drugs accurately.

### Preparation of standard stock solution of METFORMIN

10 mg of METFORMIN was weighed in to 100ml volumetric flask and dissolved in Methanol and then dilute up to the mark with methanol and prepare  $10\mu\text{g/ml}$  of solution by diluting 1ml to 10ml with methanol, wavelength is found to be 232nm.


### Preparation of standard stock solution of CANAGLIFLOZIN

10 mg of CANAGLIFLOZIN was weighed and transferred in to 100ml volumetric flask and dissolved in methanol and then make up to the mark with methanol and prepare  $10\mu\text{g/ml}$  of solution by diluting 1ml to 10ml with methanol, wavelength is found to be 290nm.

## RESULTS

The wavelength of maximum absorption ( $\lambda_{max}$ ) of the drug,  $10\mu\text{g/ml}$  solution of the drugs in methanol were scanned using UV-Visible spectrophotometer within the wavelength region of 200–400 nm against methanol as blank. The resulting spectra and the absorption curve shows the isobestic point was found to be 251 nm for the combination.

## Isobestic point of METFORMIN AND CANAGLIFLOZIN




## METHOD DEVELOPMENT OF METFORMIN AND CANAGLIFLOZIN

### Trial- 1

#### Preparation of mixed standard solution

Weigh accurately 10 mg of METFORMIN and CANAGLIFLOZIN in 100 ml of volumetric flask

and dissolve in 10ml of mobile phase and make up the volume with mobile phase. From above stock solution 10 $\mu$ g/ml of METFORMIN and CANAGLIFLOZIN is prepared by diluting 1ml to 10ml with mobile phase. This solution is used for recording chromatogram.



**Fig. 8.3.5:** Chromatogram of CANAGLIFLOZIN and METFORMIN by using mobile phase

### Observation

- Peak Asymmetry factor for CANAGLIFLOZIN and METFORMIN meet the system suitability requirements.

- The run time is very correct.
- Theoretical plates were more than 2000.
- Hence it is taken for optimization.

**Table 1: Optimized chromatographic conditions**

Mobile phase	k <sub>2</sub> hpo <sub>4</sub> +NaHPO <sub>4</sub> Buffer:ACN:MEOH (40:40:20)
Ph	5.0
Column	Inertsil ODS 3V column,C18(150x4.6 ID) 5µm
Flow rate	1.0 ml/min
Column temperature	Room temperature(20-25°C)
Sample temperature	Room temperature(20-25°C)
Wavelength	251
Injection volume	20 µl
Run time	6 min
Retention time	About 2.237 min for METFORMIN and 3.263 min for CANAGLIFLOZIN.

## ASSAY

### Preparation of samples for Assay

#### Preparation of mixed standard solution

Weigh accurately 10mg of METFORMIN and 5 mg of CANAGLIFLOZIN in 100 ml of volumetric flask and dissolve in 10ml of mobile phase and make up the volume with mobile phase. From above stock solution 10 µg/ml of METFORMIN and CANAGLIFLOZIN is prepared by diluting 1ml to 10ml with mobile phase. This solution is used for recording chromatogram.

#### Tablet sample

10 tablets (each tablet contains CANAGLIFLOZIN-50 mg METFORMIN-500 mg) were weighed and taken into a mortar and crushed to fine powder and uniformly mixed. Tablet stock

solutions of CANAGLIFLOZIN and METFORMIN (µg/ml) were prepared by dissolving weight equivalent to 10 mg of CANAGLIFLOZIN and METFORMIN and dissolved in sufficient mobile phase. After that filtered the solution using 0.45-micron syringe filter and Sonicated for 5 min and dilute to 10ml with mobile phase. Further dilutions are prepared in 5 replicates of 10µg/ml of CANAGLIFLOZIN and METFORMIN was made by adding 1 ml of stock solution to 10 ml of mobile phase.

#### Calculation

The amount of CANAGLIFLOZIN and METFORMIN present in the formulation by using the formula given below, and results shown in above table:

$$\% \text{ Assay} = \frac{AT}{AS} \times \frac{WS}{DS} \times \frac{DT}{WT} \times \frac{P}{100} \times \frac{AW}{LC} \times 100$$

Where,


AS: Average peak area due to standard preparation

AT: Peak area due to assay preparation


WS: Weight of CANAGLIFLOZIN /METFORMIN in mg

WT: Weight of sample in assay preparation


DT: Dilution of assay preparation



**Fig:** Chromatogram of Assay standard preparation-1



**Fig:** Chromatogram of Assay standard preparation-2



**Fig:** Chromatogram of Assay standard preparation-3




Fig: Chromatogram of Assay standard preparation-4




Fig: Chromatogram of Assay standard preparation-5







Fig: Chromatogram of Assay sample preparation-1




**Fig:** Chromatogram of Assay sample preparation-2



**Fig:** Chromatogram of Assay sample preparation-3



**Fig:** Chromatogram of Assay sample preparation-4



**Fig:** Chromatogram of Assay sample preparation-5

METFORMIN	CANAGLIFLOZIN			
	Standard Area	Sample Area	Standard Area	Sample Area
<b>Injection-1</b>	1487.1	1487.447	410.632	409.505
<b>Injection-2</b>	1483.265	1488.48	404.609	409.228
<b>Injection-3</b>	1488.429	1483.804	409.31	408.099
<b>Injection-4</b>	1489.131	1488.429	411.211	409.31
<b>Injection-5</b>	1483.538	1487.1	407.001	405.125
<b>Average Area</b>	1486.293	1487.052	408.5526	408.2534
<b>Assay(%purity)</b>	100.051094		99.9267659	

### Observation

The amount of CANAGLIFLOZIN and METFORMIN present in the taken dosage form was found to be 100.05% and 99.92% respectively.

### VALIDATION

#### Specificity by Direct comparison method

There is no interference of mobile phase, solvent and placebo with the analyte peak and also the peak purity of analyte peak which indicate that the method is specific for the analysis of analytes in their dosage form.

#### Preparation of mixed standard solution


Weigh accurately 10mg of METFORMIN and 10 mg of CANAGLIFLOZIN in 100 ml of volumetric flask and dissolve in 10ml of mobile phase and make up the volume with mobile phase. From above stock solution 10 $\mu$ g/ml of METFORMIN and

CANAGLIFLOZIN is prepared by diluting 1ml to 10ml with mobile phase. This solution is used for recording chromatogram.


#### Tablet sample

10 tablets (each tablet contains CANAGLIFLOZIN– 50mg METFORMIN -500 mg) were weighed and taken into a mortar and crushed to fine powder and uniformly mixed. Tablet stock solutions of CANAGLIFLOZIN and METFORMIN ( $\mu$ g/ml) were prepared by dissolving weight equivalent to 10 mg of CANAGLIFLOZIN and 20 mg of METFORMIN and dissolved in sufficient mobile phase. After that filtered the solution using 0.45-micron syringe filter and Sonicated for 5 min and dilute to 10ml with mobile phase. Further dilutions are prepared in 5 replicates of 10 $\mu$ g/ml of CANAGLIFLOZIN and METFORMIN was made by adding 1 ml of stock solution to 10 ml of mobile phase.






**Fig:** Blank chromatogram for specificity by using mobile phase



**Fig:** Chromatogram for specificity of CANAGLIFLOZIN and METFORMIN sample



**Fig:** Chromatogram for Specificity of CANAGLIFLOZIN and METFORMIN standard

**Observation**

It is observed from the above data; diluents or excipients peaks are not interfering with the CANAGLIFLOZIN and METFORMIN peaks.


**Linearity and range**

**Preparation of standard stock solution**


Standard stock solutions of METFORMIN and CANAGLIFLOZIN (microgram/ml) were prepared by dissolving 10 mg of METFORMIN and CANAGLIFLOZIN dissolved in sufficient mobile phase and dilute to 100 ml with mobile phase. Further dilutions were given in the table

**Table:** Linearity Preparations


Preparations	Volume from standard stock transferred in ml	Volume made up in ml (with mobile phase)	Concentration of solution( $\mu\text{g/ml}$ )	
			METFORMIN	CANAGLIFLOZIN
<b>Preparation 1</b>	0.75	0.24	10	75
<b>Preparation 2</b>	1	0.32	10	100
<b>Preparation 3</b>	1.25	0.4	10	125
<b>Preparation 4</b>	1.50	0.48	10	150
<b>Preparation 5</b>	1.75	0.56	10	175




**Fig:** Chromatogram of CANAGLIFLOZIN and METFORMIN preparation-1




**Fig:** Chromatogram of CANAGLIFLOZIN and METFORMIN preparation-2



**Fig:** Chromatogram of CANAGLIFLOZIN and METFORMIN preparation-3



**Fig:** Chromatogram of CANAGLIFLOZIN and METFORMIN preparation-4



**Fig:** Chromatogram of CANAGLIFLOZIN and METFORMIN for preparation-5


### Linearity of METFORMIN

S.No.	Conc.(µg/ml )	Area
1	75	927.341
2	100	1235.752
3	125	1477.14
4	150	1722.958
5	175	2041.082


### Linearity of CANAGLIFLOZIN

S.No.	Conc.(µg/ml )	Area
1	24	281.784
2	32	372.204
3	40	438.317
4	48	520.564
5	56	627.846

### Linearity graph of METFORMIN



### Linearity graph of CANAGLIFLOZIN



### Acceptance criteria

The relationship between the concentration of METFORMIN and CANAGLIFLOZIN and area of METFORMIN and CANAGLIFLOZIN should be


linear in the specified range and the correlation should not be less than 0.99.

### Observation


The correlation coefficient for linear curve obtained between concentration vs. Area for standard preparations of METFORMIN and CANAGLIFLOZIN is 0.999 and 0.996. The relationship between the concentration of METFORMIN and CANAGLIFLOZIN and area of METFORMIN and CANAGLIFLOZIN is linear in the range examined since all points lie in a straight line and the correlation coefficient is well within limits.

### Accuracy

Accuracy of the method was determined by Recovery studies. To the formulation (pre analyzed sample), the reference standards of the drugs were added at the level of 50%, 100%, 150%. The recovery studies were carried out three times and the percentage recovery and percentage mean recovery were calculated for drug is shown in table. To check the accuracy of the method, recovery studies were carried out by addition of standard drug solution to pre-analyzed sample solution at three different levels 50%, 100%, 150%



**Fig:** Chromatogram of 50% recovery (injection 1)



**Fig:** Chromatogram of 100% recovery (injection 2)




Fig: Chromatogram of 150% recovery (injection 3)




Fig: Chromatogram of 50% recovery (injection 1)




Fig: Chromatogram of 100% recovery (injection 2)




Fig: Chromatogram of 150% recovery (injection 3)




Fig: Chromatogram of 50% recovery (injection 1)




Fig: Chromatogram of 100% recovery (injection 2)




Fig: Chromatogram of 150% recovery (injection 3)

Recovery level	Accuracy METFORMIN			Average % Recovery
	Amount taken(mcg/ml)	Area	%Recovery	
100	75	410.811		102.83
	75	411.371	197.6411668	
	75	411.224	197.9105828	
			197.8398611	
120	125	522.084	125.5873028	99.36
	125	527.341	126.8518741	
	125	529.651	127.4075446	
			127.4075446	
140	175	621.351	89.67958741	88.80105072
	175	621.351	89.67958741	
	175	615.264	89.67958741	
			89.67958741	

**Acceptance criteria:** The % recovery of CANAGLIFLOZIN and METFORMIN should lie between 98% and 110%.

Table: Recovery results for CANAGLIFLOZIN

Recovery level	Accuracy CANAGLIFLOZIN			Average % Recovery
	Amount taken(mcg/ml)	Area	%Recovery	
100	24	410.811		102.83
	24	411.371	197.6411668	
	24	411.224	197.9105828	
			197.8398611	
120	40	522.084	125.5873028	99.36
	40	527.341	126.8518741	



	<b>40</b>	<b>529.651</b>	<b>99.103</b>
		<b>127.4075446</b>	
140	<b>56</b>	<b>621.351</b>	
	<b>56</b>	<b>621.351</b>	<b>89.67958741</b>
	<b>56</b>	<b>615.264</b>	<b>89.67958741</b>
		<b>88.80105072</b>	

**Observation**

The percentage mean recovery of METFORMIN and CANAGLIFLOZIN is 99.19 % and 99.89 % respectively.

**Acceptance criteria**


The % Relative standard deviation of Assay preparations of METFORMIN AND CANAGLIFLOZIN should be not more than 2.0%.

**Precision**


**Method precision**

**Method precision**


Prepared sample preparations of METFORMIN AND CANAGLIFLOZIN as per test method and injected 6 times in to the column.




**Fig: Chromatogram of precision injection 1**




**Fig: Chromatogram of precision injection 2**



**Fig:** Chromatogram of precision injection 3



**Fig:** Chromatogram of precision injection 4



**Fig:** Chromatogram of precision injection 5




Fig: Chromatogram of precision injection 6

Table 9.5.7: Results for Method precision of METFORMIN and CANAGLIFLOZIN

METFORMIN			CANAGLIFLOZIN		
S.No.	Rt	Area	S.No.	Rt	Area
1	3.145	978370.000	1	6.211	340457
2	3.165	962064.000	2	6.224	341907
3	3.151	967422.000	3	6.212	339323.000
4	3.148	955774.000	4	6.194	339473.000
5	3.126	951906.000	5	6.168	339074
6	3.116	962532.000	6	6.170	340503.000
<b>avg</b>	3.1418	963011.333	<b>avg</b>	6.197	340122.833
<b>stdev</b>	0.0178	9297.067	<b>stdev</b>	0.023	1058.443
<b>%RSD</b>	<b>0.57</b>	<b>0.97</b>	<b>%RSD</b>	<b>0.38</b>	<b>0.31</b>

**Robustness**

**Chromatographic conditions variation**

To demonstrate the robustness of the method, prepared solution as per test method and injected at different variable conditions like using different conditions like flow rate and wavelength. System

suitability parameters were compared with that of method precision.

**Acceptance criteria**

The system suitability should pass as per the test method at variable conditions.







Fig: Chromatogram of CANAGLIFLOZIN and METFORMIN Robustness (0.8 ml/min)



**Fig:** Chromatogram of CANAGLIFLOZIN and METFORMIN for Robustness (1.2 ml/min)



**Fig:** Chromatogram of CANAGLIFLOZIN and METFORMIN for Robustness (249nm)



**Fig:** Chromatogram of CANAGLIFLOZIN and METFORMIN for Robustness (253nm)

**Result of Robustness study**

Parameter	METFORMIN		CANAGLIFLOZIN	
	Retention time(min)	Tailing factor	Retention time(min)	Tailing factor
<b>Flow Rate</b>				
<b>0.8 ml/min</b>	2.562	1.679	5.059	1.263
<b>1.2 ml/min</b>	2.148	1.678	4.235	1.264

Wavelength				
249nm	2.566	1.687	5.052	1.262
253nm	2.570	1.686	5.065	1.265

**Observation**


From the observation it was found that the system suitability parameters were within limit at all variable conditions.

**Ruggedness**


The ruggedness of the method was studied by the determining the analyst to analyst variation by performing the Assay by two different analysts

**Acceptance criteria**


The % Relative standard deviation of Assay values between two analysts should be not more than 2.0%.




**Fig:** Chromatogram of Analyst 01 standard preparation



**Fig:** Chromatogram of Analyst 01 sample preparation



**Fig:** Chromatogram of Analyst 02 standard preparation



**Fig:** Chromatogram of Analyst 02 sample preparation

### Results for Ruggedness

METFORMIN	%Assay	CANAGLIFLOZIN	%Assay
<b>Analyst 01</b>	100.5	<b>Analyst 01</b>	98.9
<b>Anaylst 02</b>	99.5	<b>Anaylst 02</b>	100.6

### BIBLIOGRAPHY

- [1]. Li Q, Cheng T, Wang Y, Bryant SH. "PubChem as a public resource for drug discovery.". *Drug Discov Today* **15** (23-24), 2010, 1052–7. doi:10.1016/j.drudis.2010.10.003. PMID 20970519. edit
- [2]. Evan E. Bolton, Yanli Wang, Paul A. Thiessen, Stephen H. Bryant. "Chapter 12 PubChem: Integrated Platform of Small Molecules and Biological Activities". *Annual Reports in Computational Chemistry* **4**, 2008, 217–241. doi:10.1016/S1574-1400(08)00012-1.
- [3]. Hettne KM, Williams AJ, van Mulligen EM, Kleinjans J, Tkachenko V, Kors JA.. "Automatic vs. manual curation of a multi-source chemical dictionary: the impact on text mining". *J Cheminform* **2** (1), 2010, 3. doi:10.1186/1758-2946-2-3. PMID 20331846. edit
- [4]. Joanne Wixon, Douglas Kell. "Website Review: The Kyoto Encyclopedia of Genes and Genomes — KEGG". *Yeast* **17** (1), 2000, 48–55. doi:10.1002/(SICI)1097-0061(200004)17:1<48::AID-YEA2>3.0.CO;2-H.
- [5]. Gaulton A, Bellis LJ, Bento AP, Chambers J, Davies M, Hersey A, Light Y, McGlinchey S, Michalovich D, Al-Lazikani B, Overington JP. "ChEMBL: a large-scale bioactivity database for drug discovery". *Nucleic Acids Res* **40** (Database issue): 2012, D1100–7. doi:10.1093/nar/gkr777. PMID 21948594. edit
- [6]. Witters LA: The blooming of the French lilac. *J Clin Invest.* 2001 Oct;108(8):1105-7. PMID 11602616

- [7]. UNGAR G, FREEDMAN L, SHAPIRO SL: Pharmacological studies of a new oral hypoglycemic drug. *Proc Soc Exp Biol Med.* 95(1), 1957, 190-2. PMID 13432032
- [8]. Lord JM, Flight IH, Norman RJ: Metformin in polycystic ovary syndrome: systematic review and meta-analysis. *BMJ.* 25, 327(7421), 2003, 951-3. PMID 14576245
- [9]. Marchesini G, Brizi M, Bianchi G, Tomassetti S, Zoli M, Melchionda N: Metformin in non-alcoholic steatohepatitis. *Lancet.* 15, 358(9285), 2001, 893-4. PMID 11567710
- [10]. Nair S, Diehl AM, Wiseman M, Farr GH Jr, Perrillo RP: Metformin in the treatment of non-alcoholic steatohepatitis: a pilot open label trial. *Aliment Pharmacol Ther.* 20(1), 2004, 23-8. PMID 15225167
- [11]. Knox C, Law V, Jewison T, Liu P, Ly S, Frolkis A, Pon A, Banco K, Mak C, Neveu V, Djoumbou Y, Eisner R, Guo AC, Wishart DS. "DrugBank 3.0: a comprehensive resource for omics research on drugs". *Nucleic Acids Res.* 39 (Database issue): 2011, D1035–41. doi:10.1093/nar/gkq1126. PMC 3013709. PMID 21059682. Edit
- [12]. David S. Wishart, Craig Knox, An Chi Guo, Dean Cheng, Savita Shrivastava, Dan Tzur, Bijaya Gautam, and Murtaza Hassanali. "DrugBank: a knowledgebase for drugs, drug actions and drug targets". *Nucleic acids research* 36 (Database issue): D901–. doi:10.1093/nar/gkm958. PMC 2238889. PMID 18048412. Edit. 2008.
- [13]. Ghose, A.K., Viswanadhan V.N., and Wendoloski, J.J. "Prediction of Hydrophobic (Lipophilic) Properties of Small Organic Molecules Using Fragment Methods: An Analysis of AlogP and CLogP Methods". *J. Phys. Chem. A* **102**, 1998, 3762–3772. doi:10.1021/jp980230o.